



喜诺  
Genobio



**Rapid  
Test**

FungiXpert<sup>®</sup>

## Cryptococcal Capsular Polysaccharide Detection Kit (CLIA)

Product Code: FCrAg012-CLIA

**Obtain result in 60min**




## Overview

Cryptococcal infection can occur in several parts of the body, most commonly in the central nervous system and lungs. FungiXpert® Cryptococcal Capsular Polysaccharide Detection Kit (CLIA) is an effective product used for the quantitative detection of the Cryptococcal capsular polysaccharide in serum and cerebrospinal fluid (CSF). The assay can aid in the diagnosis of cryptococcosis in clinical. It is fully automated with FACIS to complete sample pretreatment and experimental testing, fully liberating the hands of laboratory physician and greatly improving detection accuracy.

## Product Characteristics

Basic Parameters	Cryptococcal Capsular Polysaccharide Detection Kit (CLIA)
Method	Chemiluminescence immunoassay
Sample Type	Serum/ CSF
Specification	12 Tests/Kit, 1 Test/Strip
Detection Instrument	Full-Automatic Chemiluminescence Immunoassay System
Model	FACIS-I
Detection Time	40 min
Detection Objects	<i>Cryptococcus spp.</i>
Stability	The kit is stable for 1 year at 2°C-8°C
Low Detection Limit	0.1 µg/L

## Product Advantages



**Saving time is saving life!**  
Obtain result within 60 min.



**Unique sample pretreatment system**  
using micron film with invention patent.



**All in one!**  
Integrated reagents and consumables together, more convenient.



**Fully automatically perform**  
experiment process & data analysis.

**Note:** This product can detect all the *Cryptococcus* serotypes, the corresponding serotypes are *Cryptococcus neoformans* (serotypes A&D) and *Cryptococcus gattii* (serotypes B&C), respectively.

Tel: +86(22)60978609 ext. 8605/8606/8608/8609

Fax: +86(22)60978609-804

E-mail: [marketing@era-bio.com](mailto:marketing@era-bio.com) / [customer@era-bio.com](mailto:customer@era-bio.com)

Web: [www.era-bio.com](http://www.era-bio.com)

 **Genobio Pharmaceutical Co., Ltd.**

Add: 2-1,2,3-101, No.2 Zone, Binlv Industrial Park, Binhai  
Tourism Area, Binhai New Area, Tianjin 300480, P. R. China

 **Tianjin Era Biology Technology Co., Ltd.**

Add: 1601, Building C2, TEDA MSD-C, No. 79, 1st Avenue, Tianjin  
Economic-Technological Development Area, Tianjin 300457,  
P. R. China



**Genobio is wholly-owned subsidiary of Era Biology Group**