


# eazyplex<sup>®</sup> SuperBug CRE

Carbapenemases! Clear confirmation in 15 min!



- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-detection
- **CE IVD**

Direct screening from rectal swab and urine in 20 minutes as well as confirmation from positive blood culture and culture plate in 15 minutes with eazyplex<sup>®</sup> CRE!

Due to excellent results the eazyplex<sup>®</sup> SuperBug CRE establishes its place in labs throughout Europe:


Tormo N. et al., ECCMID 2018  
Rödel J., et al., EJCMI 2017  
Cherkaou A. et al., SGM 2016  
Hinić V., et al., JMM 2015  
D'Inzeo T., et al., AMCLI 2015  
Findlay J., et al., JAC 2015  
García-Fernández S., et al., JAC 2015

The detection of the following resistances is included:

- KPC
- NDM
- VIM
- OXA-48
- CTX-M-1 group
- CTX-M-9 group
- OXA-181

See more under: [www.eazyplex.com/publications](http://www.eazyplex.com/publications)

Find out more and contact us under [support@eazyplex.com](mailto:support@eazyplex.com)



The **eazyplex<sup>®</sup>** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH  
Werkstraße 2  
83555 Gars Bahnhof  
Germany  
[www.eazyplex.com](http://www.eazyplex.com)

Kontakt/Contact

[info@eazyplex.com](mailto:info@eazyplex.com)  
Tel.: +49 (0) 8073/91693-50  
Fax: +49 (0) 8073/91693-55

