


Genie[®] HT

The new Genie High Throughput!
Your maximum of flexibility in your eazyplex[®] diagnostics!

- Stand alone device
- Up to 12 fully independent blocks
- Single retrofittable
- HIS/LIS connectivity
- Random access



Find out more and contact us under support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 80973/91693-50
Fax: +49 (0) 80973/91693-55

